

	Axis	3M Espe	Dentsply	SS White	Shofu	Dedeco	Brassler	Premier	Cosmodent	Miltex	Kerr-Hawe
Single Patient Composite Disc Kit		Sof-Lex XT			Super Snaps		EP System	Poli-Pro	Flexi-Disc		OptiDisc
JP Single-Use Diamond Composite Polishers			Pogo Enhance	Jazz							
JP Two Step Diamond Composite Polishers	PDQ2 Single Step/ DiaGloss Composite-Polisher			Jazz Supreme 1 Step							HLuster Plus Polishing System
JP Composite Strips		Sof-Lex Finishing & Polishing Strips					EP Strips	Compo Strip	Flexi-Strips	Dual Alo / Spectra-F Lightning	
JP Diamond Strips	NTI Diamond Finishing Strips									Micro Diamond	
JP Strip Saws	NTI Serrated Strips										
JP Separating Discs			Faskut Disc			Ultra Dura Very Thin	Ceramic Separating Discs			Joe Dandy	
JP Amalgam Gold Polishers	Goldstar Polishing System				Brownie and Greenie	Midgets	AmalgaPro				
JP Porcelain Polishers	CeraGlaze/NTI Porcelain Polishing			Jazz Porcelain & Metal		Hi-glaze	CeramiPro				
JP Ceramic Abrasives							Ceramic Abrasives				
JP Silicon Acrylic Polishers	NTI Acrylic Polishers					Jumbos	AcryliPro				
JP Green and White Stones	NTI Abrasives				Dura-Green Dura-White	Green Giant White Arkansas	Supra White Supra Green				

JP FINISHING & POLISHING USAGE GUIDE

	MATERIALS										
	COMPOSITE	PORCELAIN	GOLD	METALS	AMALGAM	TITANIUM	ZIRCONIA	ACRYLICS	SPRUE REMOVAL	ENAMEL	CERAMIC
POLISHING ABRASIVES											
COMPOSITE DISCS AND STRIPS	BEST										BEST
DIAMOND STRIPS	GOOD	GOOD					GOOD			GOOD	GOOD
AMALGAM/GOLD POLISHERS			BEST	BEST	BEST	BEST					
PORCELAIN POLISHERS		BEST		GOOD			BEST				BEST
SINGLE-USE COMPOSITE POLISHERS	GOOD									GOOD	
2-STEP DIAMOND COMPOSITE POLISHERS	GOOD									GOOD	
ACRYLIC POLISHERS								BEST			
SEPARATING DISCS									BEST		
FINISHING ABRASIVES											
GREEN AND WHITE STONES	BEST	BEST	BEST	BEST	BEST	BEST	BEST	GOOD	BEST	BEST	GOOD
CERAMIC ABRASIVES							BEST				BEST

Description	Item #
-------------	--------

Single Patient Composite Disc Kit 3/8"	SP-SL1981
Single Patient Composite Disc Kit 1/2"	SP-SL1982

Single Patient Composite Disc Kit helps reduce cross contamination and is conveniently packaged in individual, single patient units. Packaged 48/box, each unit contains four color coded discs and a disposable mandrel. Grits contained are one each of Coarse, Medium, Fine and Superfine.

Single-Use Diamond Composite Polishing Kit - ENH/POG-KIT

10 ea. cup, point, disc in 2 grits. Total - 60 pieces. This disposable polishing system is ideal for the finishing and polishing of all composites, especially the industry's latest generation of composites. Polishers contain diamond abrasives in two different grits.

Coarse Disc - ENH-DISC	Fine Disc - POG-DISC
Coarse Cup - ENH-CUP	Fine Cup - POG-CUP
Coarse Point - ENH-POINT	Fine Point - POG-POINT

Two-Step Diamond Composite Polishing Kit - COMP-KIT

1 ea. sm point, point, cup, knife edge in 2 grits. Total - 8 pieces. Two-step diamond polishing system for finishing, smoothing and high-luster polishing of the latest generation of composites. Diamond abrasives are embedded in the silicone polishers in two different grits.

Coarse Sm. Point - COMP-SPC	Fine Sm. Point - COMP-SPF
Coarse Point - COMP-PC	Fine Point - COMP-PF
Coarse Cup - COMP-CC	Fine Cup - COMP-CF
Coarse Knife Edge - COMP-KEC	Fine Knife Edge - COMP-KEF

Complete Disc/Strip Kit - SL2385-KIT

200 discs - 25 each 3/8" and 1/2" in Superfine, Fine, Medium and Coarse grits; 50 strips - 25 each Regular Coarse/Medium and Regular Fine/Superfine; 2 RA mandrels. Total - 252 pieces. Easy to use and ideal for all composites. Produce surface smoothness and visual gloss for more natural and life-like results. Color-coded to simplify your polishing procedures.

3/8" Coarse - SL-1981/1958C	1/2" - Coarse - SL-1982/1958C
3/8" Medium - SL-1981/1958M	1/2" - Medium - SL-1982/1958M
3/8" Fine - SL-1981/1958F	1/2" - Fine - SL-1982/1958F
3/8" Superfine - SL-1981/1958SF	1/2" - Superfine - SL-1982/1958SF

Finishing & Polishing Strips

Regular Coarse/Medium Strips	SL1954-C/M
Narrow Coarse/Medium Strips	SL1954N-C/M
Regular Fine/Superfine Strips	SL1956F-S/F
Narrow Fine/Superfine Strips	SL1956NF-S/F
Diamond Strips Super Fine	FS4SF
Diamond Strips Fine	FS4F
Diamond Strips Medium	FS4M
Perforated Diamond Strips Super Fine	FS4PSF
Perforated Diamond Strips Fine	FS4PF
Perforated Diamond Strips Medium	FS4PM

Diamond Paste

Containing actual diamond particles for a superb, natural appearance. 2gm per pkg.

Extra Fine - 520005	Regular - 520007
Fine - 520006	Assorted - 520004
8mm Strip Saw - SAW-8MM	
RA Mandrel - MAN-1958RA	HP Mandrel - MAN-1958HP
303 Mandrel - MAN-303HP	

Separating Discs - 7/8" X .009	ST-SEPD09
Separating Discs - 7/8" X .015	VT-SEPD15
Separating Discs - 7/8" X .025 - Ruby	RUBY-DO25
Separating Discs - 7/8" X .025 - Stripper	STRIP-DO25

Amalgam/Gold Polishing Kit - AM-GOLDKIT

1 ea. mini point, point, cup, bullet in 3 grits for polishing amalgam and metal. Total 12 pcs. Three-step polishing system for removing, smoothing and polishing amalgam and gold. The silicon carbide is embedded in a medium-soft, flexible synthetic rubber matrix.

Shaping/Coarse Mini Pt-AMAL-P030C	Shaping/Coarse Cup-AMAL-C065C
Shaping/Coarse Point-AMAL-MP050C	Shaping/Coarse Bullet-AMAL-B040C
Smthing/Med Mini Pt-AMAL-MP050C	Smthing/Med Cup-AMAL-C065M
Smthing/Med Point-AMAL-MP050M	Smthing/Med Bullet-AMAL-B040M
High Shine/Fine Mini Pt-AMAL-MP030F	High Shine/Fine Cup-AMAL-C065F
High Shine/Fine Point-AMAL-MP050F	High Shine/Fine Bullet-AMAL-B040F

Description	Item #
-------------	--------

Diamond Porcelain Polishing Kit - CERA506-KIT

1 ea. sm point, point, cup, disc in 3 grits for polishing of ceramic, metal and zirconia. Total-12 pieces. Three-step diamond polishing system for removing, finishing and high-luster polishing of ceramic, metal and zirconia. Diamond abrasive is embedded in a hard- lexible, silicone polisher. The high-luster polish is achieved without paste.

Shaping/Coarse Wheel - CERA-W43C	Smthing/Med Wheel - CERA-W43M
Shaping/Coarse Sm Pt - CERA-SP33C	Smthing/Med Sm Pt - CERA-SP33M
Shaping/Coarse Point - CERA-P32C	Smthing/Med Point - CERA-P32M
Shaping/Coarse Cup - CERA-C55C	Smthing/Med Cup - CERA-C55M
Shaping/Coarse Disc - CERA-D36C	Smthing/Med Disc - CERA-D36M
Hi Shine/Fine Wheel - CERA-W43F	Hi Shine/Fine Cup - CERA-C55F
Hi Shine/Fine Sm Pt - CERA-SP33F	Hi Shine/Fine Disc - CERA-D36F
Hi Shine/Fine Point - CERA-P32F	

Standard Porcelain Polishing Kit - PORC-KIT

1 ea. sm point, point, cup, disc in 3 grits for polishing of porcelain and enamel. Total-12 pieces. Silicone polishers made specifically for polishing porcelain and enamel, with three grits for pre polishing, polishing and high shine super polishing.

Shaping/Coarse Wheel - PORC-W43C	Smthing/Med Wheel - PORC-W43M
Shaping/Coarse Sm Pt - PORC-SP33C	Smthing/Med Sm Pt - PORC-SP33M
Shaping/Coarse Point - PORC-P32C	Smthing/Med Point - PORC-P32M
Shaping/Coarse Cup - PORC-C55C	Smthing/Med Cup - PORC-C55M
Shaping/Coarse Disc - PORC-D36C	Smthing/Med Disc - PORC-D36M
Hi Shine/Fine Wheel - PORC-W43F	Hi Shine/Fine Cup - PORC-C55F
Hi Shine/Fine Sm Pt - PORC-SP33F	Hi Shine/Fine Disc - PORC-D36F
Hi Shine/Fine Point - PORC-P32F	

Ceramic Abrasives Kit - ZIRC-KIT

1 ea. disc, cone in 2 grits. Total - 4 pieces. High performance synthetically bound diamond abrasives for removing metal, specially designed for zirconia and alumina restorations without generating heat. Coarse abrasives for gross reduction, medium abrasives for finishing. For polishing, use Porcelain Polishers.

Coarse Disc - ZIRC-DC	Coarse Cone - ZIRC-CC
Med Disc - ZIRC-DM	Med Cone - ZIRC-CM

Silicone Acrylic Polishing Kit - ACR-KIT

1 ea. sm point, large point, jumbo in 3 grits. Total-12 pieces. Three-step polishing system for removing, finishing and high luster polishing of acrylics. Silicon carbide embedded in a medium-soft, flexible synthetic rubber matrix.

Shapin Crse Sm Pt - ACR-SP9580C	Shapin Crse Lg Pt - ACR-LP674C
Shapin Crse Point - ACR-P679C	Shapin Crse Jumbo - ACR-J672C
Smthing Med Sm Pt - ACR-SP9580M	Smthing Med Lg Pt - ACR-LP664M
Smthing Med Point - ACR-P669M	Smthing Med Jumbo - ACR-J662M
Hi Shine Fine Sm Pt- ACR-SP9580F	Hi Shine Fine Lg Pt - ACR-LP654F
Hi Shine Fine Point - ACR-P659F	Hi Shine Fine Jumbo - ACR-J652F

Green & White Stones

Green Stones contain silicon carbide and are excellent for fast contouring and finishing of porcelain, composites, and gold and silver alloys. White Stones, made of aluminum oxide, are ideal for fine contouring and finishing of enamel, composites, and porcelains.

White Stone WH-6 - W-WH6HP-12	White Stone IC-4 - W-IC4HP-12
White Stone TC-4 - W-TC4HP-12	White Stone CN-1 - W-CN1FG-12
White Stone FL-3 - W-FL3HP-12	White Stone FL-2 - W-FL2FG-12
Green Stone IC-3 - G-IC3HP-12	Green Stone WH-6 - G-WH6HP-12
Green Stone IC-5 - G-IC5HP-12	Green Stone CY-3 - G-CY3HP-12
Green Stone IC-7 - G-IC7HP-12	Green Stone TC-2 - G-TC2HP-12
Green Stone IC-9 - G-IC9HP-12	Green Stone TC-4 - G-TC4HP-12
Green Stone KN-9 - G-KN7HP-12	Green Stone CN-1 - G-CN1FG-12
Green Stone PC-2 - G-PC2HP-12	Green Stone FL-2 - G-FL2FG-12
Green Stone WH-3 - G-WH3HP-12	