

High Speed Handpieces

Gemini Series

Comparable design and performance to the market leader, but at a much better value

- 3-port misting water spray for excellent bur cooling and more efficient cutting
- Push button autochuck for quick bur changes
- Ceramic bearing construction for quiet and long lasting performance
- High power torque of 30 watts, operating at 400,000 rpm
- Stainless steel head and anti-slip surface
- Lightweight and ergonomic design weighing only 45 grams
- Safety features include safety quick stop turbine and anti-retraction valve to avoid cross contamination
- Fully autoclavable



Bella Torque Fixed Back High Speed

G-KBTFBPB

- Bella Torque's smaller head provides easier access and quiet operation (only 61 dB) (13mm height x 11mm diameter)
- available in 2 hole and 4 hole configuration

Super Torque Fixed Back High Speed

G-KSTFBPB

- Super Torque provides industry-leading torque for fast, efficient preparations and only 60 dB sound (14mm height x 12mm diameter)
- available in 2 hole and 4 hole configuration