High Speed Handpieces

Gemini Series

- Gemini handpieces offer exactly the same design and performance as the originals for a fraction of the price!
- Gemini handpieces fit the coupling systems of the originals and offer complete interchangeability with the original handpieces, turbines, and component parts



- · Provides smooth and controlled cutting with high torque
- · Durable, even after repeated sterilization
- · Mid-sized head for superior access with high torque
- · Bright twin beam fiber optics eliminates shadows
- · Fine mist water spray
- ·Lightweight

	PART #	CROSS REFERENCE #
Non-Fiber Optic		
Standard Tool Type	G-MWT	sc#7720733, 7721898, mw#750044
Push Button	G-MWTPB	sc#7726063,7724247, mw#780044,790044
Fiber Optic		
Standard Tool Type	G-MWTFO	sc#7729466, mw#750045
Push Button	G-MWTFOPB	sc#7720887, 7727345, mw#780045, 790045
FUSITEUTION	G-MWIFOPB	



High Speed Handpieces

Gemini Series

- Gemini handpieces offer exactly the same design and performance as the originals for a fraction of the price!
- Gemini handpieces fit the coupling systems of the originals and offer complete interchangeability with the original handpieces, turbines, and component parts

KAVO TYPE

- · Ceramic bearings for smooth, quiet operation
- · Durable glass cellular optics providing excellent brightness
- · Multiflex LUX connection
- · Push button auto chuck
- · 3 port water spray for excellent bur cooling and more efficient cutting
- · Super torque provides industry-leading torque for fast, efficient preparations
- · Bella Torque's smaller head provides easier access and quiet operation



Push Button Only Bella-Torque Style (smaller head, easier access)	PART #	CROSS REFERENCE #
Non-Fiber Optic	G-KBTPB	Kavo 632/636,642,643CP sc#6284588
Fiber Optic	G-KBTFOPB	Kavo 632/642/635/4500
Super-Torque Style (largest head, strongest torque)		sc#6282596,br#5012737u0
Non-Fiber Optic	G-KSTPB	Kavo 625CD,646 sc#6284885
Fiber Optic	G-KSTFOPB	Kavo 625/640/647/649/659/6000/6500 sc#6280504,br#5012738u0

Kavo Multiflex-type Couplers

- · Enables quick instrument changes in the operatory
- · 360° full rotation for improved ergonomics and safety
- · Built-in anti-retraction valve prevents the backflow of potentially contaminated water in the tubing
- · Quick and convenient connection to remove handpieces for cleaning and sterilization
- · LED option provides unmatched illumination, and light source lasts 10 times longer than halogen bulbs



Kavo Multiflex-type Couplers	PART #	CROSS REFERENCE #
Non-Optic 4 Hole	W-KQDS	457 sc#6287944
Fiber Optic 5 Hole	W-KFOQDS	1390 sc#6289674,br#5014226u0
Fiber Optic 6 Pin	W-K6PQDS	465LRN, 460LE sc#6287519,6282022
Fiber Optic 6 Pin LED	W-K6PQDS-LED	460LED,465LED sc#6280337,6280339



High Speed Handpieces

Gemini Series

- Gemini handpieces offer exactly the same design and performance as the originals for a fraction of the price!
- Gemini handpieces fit the coupling systems of the originals and offer complete interchangeability with the original handpieces, turbines, and component parts



- · 360° swivel to reduce fatigue and for convenient quick-click connections
- · Mid-size head provides excellent balance of high torque and good intra-oral access

Push Button Only	PART #	CROSS REFERENCE #	
Non-Fiber Optic	G-VSPB	N75 br#500980100	
Fiber Optic	G-VSFOPB	Viper 360,NL85 br#5009781u0,5012734u0, 5012733u0	

	NSK-type Couplers	PART #	CROSS REFERENCE #
4 HOLE	Non-Optic 4 Hole	W-NQDS	NL11 br#50098200,5009809u0
5 HOLE	Fiber Optic 5 Hole	W-NFOQDS	N10 br#5009875u0,5009802u0
6-PIN	Fiber Optic 6 Pin	W-N6PQDS	NL10 br#5009820u0,5009820u0

