



Johnson-Promident

Clear Blue Digi 2.0 LED Curing Light

*One of the earliest
LED curing lights
and as popular
as ever!*

- High power 1,200 mW/cm² for fast cures
- Emits light in the most relevant spectrum of 430-480 nm wavelength to efficiently solidify any brand of resins
- Three working modes and ability to set exposure times for any clinical protocol
- Super-compact, lightweight, and cordless for the ultimate in convenience
- Cool and silent; no fan necessary as with halogen lamps
- Li-Ion battery with super capacity and no memory effect
- Digital display of working time
- Can also be operated with the cord attached for back-up in case of low battery charge
- One year warranty

