

Disposable, Evacuation and Infection Control Products CATALOG



Nitracare





JP Solutions and its family of brands – **Weltex, Ecotex, Medgluv, Neuskin, Nitracare, and Dazzle** – are proud to present our newest product catalog.

We strive to offer premium quality products at prices that will allow you to provide excellent value to your customers while maintaining strong margins.

Why Ecotex is Unique:



We are deeply concerned about our environment and realize that disposable products have a negative impact on the planet. As a supplier of these products, we also understand the critical functions they play in a safe and productive clinical environment. In an effort to balance the needs of the planet and the practitioner, we have begun researching and developing **Ecotex** brand products that provide the form and function that your customers require but integrate the responsible application of the 3 Rs (Reduce, Reuse, Recycle).

Many changes are meaningful yet incremental as we search for new technologies and materials to continue to provide more environmentally friendly products. **Ecotex** products have innovative packaging designs or methods to **REDUCE** the use of packaging materials and **REDUCE** oil consumption in transportation. **Ecotex's** current dental sponge offering can save 500 tons of paper for the industry.

Ecotex products and packaging make use of natural, biodegradable and **RECYCLED** materials where possible to save trees and reduce landfill impact.

The name **Ecotex** is derived from Ecology and Economy, because these positive environmental efforts actually cost less. As we strive to bring more **Ecotex** products to market and reduce the impact of those currently available we will donate 5% of all profits from **Ecotex** sales to non-profit conservancy organizations.

TABLE OF CONTENTS

Disposable Products

Cotton Rolls.....	3
Patient Towels/Bibs.....	3
Non-Woven Sponges.....	3,4
Disposable Prophyl Angles.....	4
Prophyl Paste.....	4

Evacuation Products

Saliva Ejectors.....	5
Evacuation Tips.....	5
Surgical Aspirator Tips.....	5
Disposable Air/Water Syringe Tips.....	5
Evacuation System Cleaner.....	5
Air/Water Syringe.....	5
Syringe Tips.....	5
Quick Disconnect Vacuum Valve.....	5
Quick Disconnect Saliva Ejector Valve.....	5

Infection Control Products

Powder-Free Latex Gloves.....	6
Powder-Free Nitrile Gloves.....	6
Powder-Free Vinyl Gloves.....	6
Sterilization Pouches.....	6
Hard Surface Disinfectant.....	6
Ear Loop Face Masks.....	7
Surgical Tie-On Face Masks.....	7
Barrier Film.....	7
Syringe Barrier Sleeves.....	7
Tray Sleeves.....	7
Headrest Covers.....	7
Light Handle Sleeves (T-Style).....	7
Chair Sleeves.....	7

Disposable Products

WELTEX

WELTEX Cotton Rolls

- Made from soft, super-absorbent 100% cotton
- Lint-free and smooth
- Maintain shape when wet and will not adhere to mucous membranes
- Pure white
- Non-sterile



#180411 Medium #2 (1½"x3")
2,000 pc/box
12 boxes/case

WELTEX Non-Woven Sponges

- Premium 50 gram weight is most absorbent non-woven sponge available
- Made from a 4-ply rayon/polyester blend which provides for maximum absorbency and minimizes adhesion to wounds
- Significantly greater absorbency for the weight than cotton
- Soft, strong, smooth
- Virtually lint-free so no fibers are left in the patient's mouth
- Larger sizes ideal for surface and instrument cleaning
- Non-sterile



200 pc/sleeve
25 sleeves/case
#180103 Premium 2x2/4 ply/50g

200 pc/sleeve
10 sleeves/case
#180303 Premium 4x4/4 ply/50g

Call for custom quotes on 3x3/4 ply/50 g

WELTEX Patient Bibs/Towels

- The unique technique used to bond the layers of tissue and poly securely together eliminates separation of the layers
- Reinforced edges provide added strength and durability
- Heavy textured 2-ply tissue and 1-ply poly for absorption and patient protection
- Available in 5 popular colors



13"x19" 500 pc/case
#430100 Blue
#430200 Lavender
#430300 Dusty Rose
#430400 White
#430500 Green

Disposable Products

ECOTEX

ECOTEX Non-Woven Sponges Premium/50gram

- Made from a 4-ply rayon/ polyester blend which provides for maximum absorbency and minimizes adhesion to wounds
- Premium 50 g weight is the most absorbent non-woven available
- Soft, strong, smooth
- Virtually lint-free so no fibers are left in the patient's mouth
- Larger sizes ideal for surface and instrument cleaning
- Non-sterile
- Manufactured with the environment in mind, including the use of recycled packaging and biodegradable inks. Special compressed packaging produces a 10% reduction in packaging material and a reduction in fuel consumption. Significant shipping savings, as 15% more product will fit on a pallet.



200 pc/sleeve
25 sleeves/case
#22100 Premium 2x2/4 ply/50g
Call for custom quote on 3x3/4 ply/50g
and 4x4/4ply/50g

ECOTEX Non-Woven Sponges Standard/35 gram

- Made from a 4-ply rayon/ polyester blend which provides for maximum absorbency and minimizes adhesion to wounds
- Soft, strong, smooth
- Virtually lint-free so no fibers are left in the patient's mouth
- Larger sizes ideal for surface and instrument cleaning
- Non-sterile
- Manufactured with the environment in mind, including the use of recycled packaging and biodegradable inks. Special compressed packaging produces a 15% reduction in packaging material and a 10% reduction in fuel consumption. Significant shipping savings, as 15% more product will fit on a pallet.



#22300 Standard 2x2/4 ply/35g
 200 pc/sleeve 25 sleeves/case
#22400 Standard 3x3/4 ply/35g
 200 pc/sleeve 20 sleeves/case
#22500 Standard 4x4/4 ply/35g
 200 pc/sleeve 10 sleeves/case

ECOTEX Disposable Prophy Angles

- Ultra-strong unit construction to reduce vibrations
- Latex-free webbed and ribbed cup
- Bevel gear design provides smooth, consistent operation
- Minimizes splatter and vibration
- Plastic, individually sealed bags
- With firm cup with individual stripes; do not exceed 3,000 RPM
- Adapt to all types of low speed handpieces (Doriot or E-type)
- Manufactured with the environment in mind. Package uses recycled materials and biodegradable inks and has a unique space-saver design that reduces the use of packaging materials by about 40%, significantly reducing shipping costs and fuel used in transportation.



100 pc/box **#23103 Firm**
6 boxes/case **#23101 Soft**

DAZZLE

DAZZLE Prophy Paste

- Non-splattering, gluten-free prophy paste
- Rinses completely and easily
- Contains 1.23% fluoride
- Great tasting

200 pc/box
12 boxes/case



#610100 Extra Coarse Mint
#610110 Coarse Mint
#610120 Medium Mint
#610130 Fine Mint
#610210 Coarse Bubblegum
#610220 Medium Bubblegum
#610230 Fine Bubblegum
#610310 Coarse Raspberry
#610320 Medium Raspberry
#610910 Coarse Assorted
#610920 Medium Assorted



Evacuation Products

WELTEX / ECOTEX

WELTEX Saliva Ejectors

- Fit all standard saliva ejector hose ends/6" length
- Bends to the right curve for maximum patient comfort
- High flow, smooth, flexible bonded tips won't come off
- Slotted bullet design of tip prevents suction of tissue and increases patient safety and comfort
- No memory; tip holds bend

100 pc/bag

10 bags/case

#850233 Clear with Blue Tip

#850231 White with White Tip




WELTEX Evacuation Tips

- Fit all standard evacuator hose ends
- One end vented, one end non-vented (dual use)
- Soft rigid tubing allows for greater patient comfort yet is strong enough for cheek retraction
- Smooth edges for maximum patient comfort
- Thoroughly polished surface for smooth and quiet operation
- Straight handle allows for increased control
- Convenient single-use disposable tips

100 pc/bag

10 bags/case

#850241 Vented/Non-Vented



WELTEX Surgical Aspirator Tips

- Made with smooth edges for maximum patient comfort
- Surgical aspirator is molded at a 30 degree angle and fits 11mm suction holders
- Autoclavable up to 135 C or may be disposed after single use

25 pc/bag

20 bags/case

#850251 1/4" Green

#850253 1/8" White



WELTEX Disposable Air/Water Syringe Tips


- Fit all standard air/water syringes without adapter
- Interchangeable with metal tips
- No leakage, guaranteed dry air
- Separate air and water tubes eliminate cross-over
- Smooth edge

250 pc/bag

6 bags/case

#850261 White

#850265 Assorted Colors



ECOTEX Evacuation System Cleaner

- Bacterial-based cleaner provides residual over-night debris removal from cuspidor drains and suction systems
- Works under all conditions, even in cold water
- Neutral pH, safe for amalgam traps
- Residual bacterial action; bacteria collect in bends and turns where they continue to liquefy debris
- Safe for equipment, environment and staff
- Pleasant scent
- Incredibly economical; 128 treatments per gallon
- Environmentally friendly; virtually foam free and 6 natural non-pathogenics work harmoniously to clean drains and suction lines

64 oz. bottle

4 bottles/case

#25301



Air/Water Syringe

- Well balanced design, and smooth styling for comfortable use
- Easy release for a speedy exchange of tips
- Fully autoclavable tips



SYR Air/Water Syringe

SYR-TIP Replacement Tip

Universal Quick Disconnect Vacuum Valve

- Single piece anodized aluminum body
- Easy removal lever and spool valve assembly
- Lightweight 360° quick disconnect swivel
- Autoclavable



HVE-QD

Deluxe Quick Disconnect Saliva Ejector Valve

- Single piece anodized aluminum body
- Easy removal lever and spool valve assembly
- Lightweight 360° quick disconnect swivel
- Autoclavable



DSE-QD



Infection Control Products

ECOTEX/NEUSKIN/NITRACARE

ECOTEX/MEDGLUV Powder-Free Latex Gloves

- Low protein for increased safety and minimum risk of protein allergies (50 mcg or less of water-extractable protein, the lowest defined content)
- Premium choice meets surgical glove 1.0 AQL standards
- Textured grip for both moist and dry environments
- Polymerized to facilitate donning
- Soft modulus for a precise feel and comfortable fit during extended wear
- Powder Free to reduce air-born particulates
- Manufactured with the environment in mind, including the use of all natural and biodegradable latex product packaging and inks. Biofuels are used for the manufacturing process. Special packaging design and packing procedure produces a 20% reduction in packaging material and a 15% reduction in fuel consumption. Significant shipping savings, as 15% more product will fit on a pallet.

100 pc/box
10 boxes/case
#211010 Extra-Small
#211020 Small
#211030 Medium
#211040 Large

NITRACARE Powder-Free Nitrile Gloves

- The protection of nitrile with the feeling of latex.
- 100% synthetic nitrile polymer is latex-free, powder-free, and odor-free.
- Unique thermal sensitive polymer compound absorbs body heat and easily conforms to user's hands to provide comfort for extended periods.
- Approved for use with chemotherapy drugs.
- Meet or exceed ASTM D6319.
- Manufacturing plants are built with clean technology for minimal impact on the environment.
- Packed 200 gloves per box to minimize storage space required.



200 pc/box
10 boxes/case

#MG-505XS Extra-Small
#MG-505S Small
#MG-505M Medium
#MG-505L Large
#MG-505XL Extra-Large

NEUSKIN Powder-Free Vinyl Gloves

- Unique advanced formulation vinyl that stretches and has soft elastic feel with a more snug fit than traditional vinyl
- Soft stretch vinyl with latex look and feel
- Special vinyl processing allows for easy donning without coatings or powders
- 100% latex free; great alternative for people with latex sensitivity



100/box
10 boxes/case
#MG-325S Small
#MG-325M Medium
#MG-325L Large
#MG-325XL Extra-Large



WELTEX Assurance Dual-Indicator Sterilization Pouches

- Internal and external indicators
- Triple-sealed seams
- Dual indicator for steam and ETO
- Medical grade paper
- Perforated fold



200/BOX
10 BOXES/CASE
#420100 2.25" X 4"
#420200 2.75" X 9"
#420300 3.5" X 5.25"
#420400 3.5" X 9"
#420500 5.25" X 10"



WELTEX Hard Surface Disinfectant

- Ready-to-use one-step general disinfectant, cleaner and deodorizer
- Effective against TB, HBV, HIV, Staph, etc.
- 100% alcohol- and phenol-free
- EPA registered
- Works on all hard surfaces
- Leaves a clean scent
- Ready-to-use



32 oz. Spray
12 bottles/case
#850125

1 Gallon Refill
4 bottles/case
#850135



Infection Control Products

WELTEX

Weltext Face Masks

- Made of a lightweight, 3-ply material that provides the ultimate in comfort with a high BFE rate
- Fluid resistant, but with excellent breathability characteristics
- Built-in adjustable nose piece maintains shape for a better fit
- Soft inner lining
- 100% latex free and fiberglass free
- Full-size pleating to completely cover nose, mouth and chin


Ear Loop

50/box
10 boxes/case
#150151 Blue
#150153 Pink




Surgical Tie-On

50/box
10 boxes/case
#150401 Blue



Weltext Barrier Film

- Universal infection control barrier with a low tack adhesive backing that is perfect for use in areas that are difficult to clean, disinfect or sterilize.
- Perforated 4"x6" sheets are easy to apply, easy to remove and leave behind no residue.



1200 pc/roll
12 rolls/case
#410200 Blue
#410210 Clear

Weltext Syringe Barrier Sleeves

- Pre-cut opening for easy application
- Economical and convenient
- Ideal for 3-way syringes, saliva ejectors, HVE valves and handpieces
- 10" x 2"

500 pc/box
15 boxes/case
#410100



NEW BARRIER PRODUCTS

Weltext barrier products are constructed with thick clear plastic material. Easy application and disposal.

Weltext Tray Sleeves

- Designed to cover the entire instrument tray

500/box
6 boxes/case
#410320 10.5" x 14"



Weltext Headrest Covers

- Designed to fit over the top of the headrest

250/box
12 boxes/case
#410370 10" x 11"
#410380 10" x 14"



Weltext Light Handle Sleeves (T-Style)

- Designed to fit most light handles

500/box
12 boxes/case
#410310 4" x 5.75"



Weltext 1/2 Chair Sleeves

- Designed to easily fit over most chairs

225/box
3 boxes/case
#410300 27.5" x 24"



Terms and Conditions

Returns: Before opening a sealed box or carton, please examine the label to ensure you have received the correct item.

Prior approval must be obtained from JP SOLUTIONS before merchandise may be returned for credit. Only standard merchandise of current stock, less than six months from original shipping date and accompanied by the original invoice will be accepted for credit. Shipment of returned merchandise must be prepaid.

Restocking fee: Merchandise approved and returned in the original unopened container and accompanied by the invoice, will be charged a 20% restocking fee.

Conditions: Purchasing of JP SOLUTIONS products from this list is open to recognized accounts and/or dealers of JP SOLUTIONS. Possession of a dealer's price list does not constitute an offer to sell. JP SOLUTIONS reserves the right to make changes in product design without any obligation to accept previously manufactured products for alteration, exchange or credit. Payment terms are net 30 days. A late payment fee of 1 1/2% will be assessed monthly on the total invoice past due.

Claims: Shortage claims must be made within 48 hours from receipt of the goods. Damaged goods must be reported to the carrier within 24 hours and damage claims must be filed by the consignee.



Disposable ■ Evacuation ■ Infection Control