




WELTEX Saliva Ejectors

- Fit all standard saliva ejector hose ends / 6" length
- Bends to the right curve for maximum patient comfort
- High flow, smooth, flexible bonded tips won't come off
- Slotted bullet design of tip prevents suction of tissue and increases patient safety and comfort
- No memory; tip holds bend
- Flexible PVC body, translucent, transparent and free of imperfections and impurities
- Internal wire made of stainless steel alloy which is covered and malleable
- PVC nozzle tip contains both lateral and central perforations and is flexible, smooth, fixed into place and has a rounded atraumatic cap
- Lumen will not be blocked when bent
- Internal diameter 5mm, External diameter 7mm



100 pc/bag

10 bags/case

#850233 Clear with Blue Tip

#850231 White with White Tip



Johnson-Promident

www.johnsonpromident.com